

Cambridge NanoTech ALD System – User Instructions

This is a shortened version that focuses mainly on the operation. For more detailed instructions, please refer to the user manual “Cambridge NanoTech Savannah 100 Atomic Layer Deposition System” available both on-line and near the ALD setup.

System description


Atomic layer deposition (ALD) is a technique that allows growth of thin films with atomic layer precision. The ALD system in PRISM clean room is designated for growth of Al_2O_3 film, which uses water and trimethylaluminum (TMA) as precursors. ALD is a self-limiting process so each cycle produces exactly a monolayer with a maximum of 1.1\AA of Al_2O_3 depending on the temperature.

Before starting: Important notes

- **TMA is pyrophoric and ignites when in contact with the air. Never remove TMA source.**
- **The chamber lid and the walls are HOT!** Use care when opening the chamber. Do not place any flammable materials on or near the ALD machine.
- **The chamber lid cannot be lifted if the chamber is cold.** The temperature of the chamber outer heater should be set to at least 80°C .
- Do not put in materials that will outgas. Do not use a substrate having carbon tape. Do not use vacuum grease on your substrate or in the chamber.
- System is intended to have the pump running, nitrogen flowing, and heaters on at all times, even when idle.
- If restarting the system after it has been cold, e.g. following a power outage or maintenance shutdown, warm up the system first and follow the out-gassing procedure completely as described in the system manual. If you are not sure, contact Helena, Joe or one of the superusers to restart the system.
- In case that either of the precursors is running out, indicated by the disappearance of pulses during the running cycles, contact Helena, Joe or one of the superusers for replacement. Do not attempt to change precursors by yourself.
- The largest wafer size supported is 4” circular.
- **Always set the center heater to 150°C and the outer heater to 150°C after you finish the growth.**

Software user interface

Buttons:



Program — Stop the Labview program

Stop Valve — Open and close the stop valve to the pump

All heaters — Turns all the heaters ON/OFF

Run cycle area:

Valve	Pulse (s)	Expo (s)	Pump(s)
0	0.1	0	5
1	0.1	0	5

Cycles	Flow (sccm)	Delay (min)
100	20	0
53	20.4	

☒ **Save Run**

C:\ALD data\ALD run Jun-11, 2005, 0546PM06sec

Runtime 00:16:59

Abort Run

Cycle matrix — Programming of the cycle recipe. **Valve** indicates the source valve to be opened. Valve 0 connects to the DI water precursor, and 1 to the TMA precursor. **Pulse** sets the time for which the source valve is opened. Valves only respond to pulsing for $t > 0.007$ sec. **Expo** gives the time for which the stop valve is closed after pulsing the source and before pumping the chamber. See exposure mode below. **Pump** sets the time between the pulses of each source, in which the chamber is pumped and vapor removed. This should be long enough so that the chamber reaches base pressure before next pulse.

Cycles — Programmed number of cycles. The current cycle number is displayed underneath.

Flow — Sets the flow rate of carrier gas (here nitrogen with ultra-low water content) as controlled by MFC.


Delay — Time delay before the start of a run. Usually set for temperature stabilization.

Save Run — Check the box to save the run data in the file described below.

Run Time — Calculated total time for the current run.


Start/Abort Run — Start or abort the current run. The arrow turns bright green as run starts and gray green as it stops.

Reactor Pressure Plot Area:



This plot tracks the reactor pressure read by the gauge installed in the pumping line assembly. Pulses can be seen on the plot. Time scale can be reset by **Reset Time**.

Heaters control Area:



Control the temperature of the heaters. Temperatures can be set by typing values into the light blue areas. The current temperature readings are indicated in the red area of each heater under the set point.

Heater settings: **Maximum** / **Recommended** (°C)

- Solenoid valve — 115 / 80
- Center heater — 400 / 250-300
- Outer heater — 200 / 150 (80 minimum)
- Tee and elbow — 180 / 150
- Pumping line — 180 / 150

Operating Instructions:

Load/Unload

1. Close the stop valve and set the flow rate to 100 sccm to vent the chamber. Only vent when system hot (chamber >60°C). Leave system vented for as short a time as possible.
2. Wait for the lid to rise slightly from the chamber. Please be patient, **it is important to wait for the chamber to reach atmospheric pressure to preserve the integrity of the O-ring.**
3. Lift the lid and place your sample in the center of the chamber. Close the lid.
4. Open the stop valve to pump down the chamber. The base pressure should be less than 1 Torr, depending on the flow rate.

Growth run —Continuous flow mode (recommended)

1. Pump down the reactor chamber by opening the stop valve after loading sample.
2. Make sure all heaters are set to the desired temperatures.
3. Enter the desired values for cycle parameters. A typical recipe is below. This recipe below will give you 300 monolayers (approx. 33 nm) of Al₂O₃.

Valve	Pulse (s)	Expo (s)	Pump (s)	Cycles	Flow (sccm)	Delay (min)
0	0.07	0	4	300	20	0
1	0.07	0	4			

4. Click **Start run** in the run cycle area. You can see the pulses in the pressure plot area.
5. After the run is over, vent the system and remove your sample. Close the lid and pump down with flow rate set to 20 sccm. Let the system idle.

Growth run —Exposure mode

Exposure mode can be used for the ultra high aspect ratio structures. A typical recipe is listed below. For both valves, **Pulse** sets to 0.2 sec; **Expo** sets to 4 sec to allow fully reaction of the current precursor with the sample surface; **Pump** sets to 10 sec to provide sufficient time to remove any un-reacted precursor and/or by products. Carrier gas is recommended and **Flow** can use 5sccm.

Degassing (only after installation of precursors or system checked for leaks)

Water cylinder — **Valve 0, Pulse = 1.3, Expo = 0, Pump = 3, Cycle = 5, Flow = 5**
TMA cylinder — **Valve 1, Pulse = 1.3, Expo = 0, Pump = 3, Cycle = 5, Flow = 5**

Keep the manual valve of TMA source still closed while degassing. Degassing is essential every time the space between manual and solenoid valves is open to air.